SAFE® Complete Care Competence

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Édito

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Comparison of two high-fat diets for inducing obesity

Clara Sanchez (sanchez@ipmc.cnrs.fr) and
Dr Carole Rovere (rovere@ipmc.cnrs.fr)
Institute of Molecular and Cellular Pharmacology (IPMC)
CNRS UMR 7275 / Université Côte d'Azur (France)

OBESITY

As a major public health problem, obesity is defined by WHO as an abnormal or excessive fat accumulation that may impair health. It's commonly used to classify overweight and obesity with a Body Mass Index (BMI) and in 2016, 39 % of adults worldwide were overweight (BMI than or equal to 25) and 13 % of adults worldwide were obese (BMI than or equal to 30). The pathology of obesity is still increasing worldwide and is linked to more deaths than underweight. Principal complications of obesity are cardiovascular diseases (mainly heart disease and stroke) which are the leading cause of death for at least 10 years. Obesity is also causing type 2 diabetes and certain types of cancer. Moreover, obesity is most commonly caused by excessive food intake and lack of physical activity making it preventable. Obesity studies are then essential to better understand the pathology and find a way to prevent or treat it and its related diseases.

MICE MODEL USED FOR NUTRITIONAL OBESITY

Mice fed healthy or harmful diets are ideal experimental models for these studies because of the way it can affect their weight and health quite similarly as human. This is why scientists working on obesity are using high-fat diets to induce the pathology in mice. We used 4 weeks old male mice from the C57BL6/J strain of Janvier breeder (St-Berthevin, France) and kept them in inverted cycle of 12 hours light and 12 hours darkness. In order to compare two different high-fat diets, we tested the U8978 version 19 High Fat Anhydrous Butter diet (SAFE® 260 HF) and the U8955 version 19 High Fat Anhydrous Butter diet (SAFE® 246 HF) for a period of the 12 weeks, compared to mice nourished on a standard SAFE® A03 diet.


Figure 1 : Effect of High Fat diets on weight gain in mice. The mice (n=10 for each group) aged 4 weeks at the beginning of the trial were weighted every week throughout the feeding with SAFE® 260 HF and SAFE® 246 HF compared to mice fed with the standard SAFE® A03 diet (up to 20 weeks). The data is \pm SEM and Kruskal-Wallis statistic test is performed at each time. * $p < 0.05$; ** $p < 0.01$; *** $p < 0.001$; **** $p < 0.0001$

C57BL6/J mice n=10 per diet	Animal's weight (g)
SAFE® A03	28,6 ± 1,6
SAFE® 246 HF U8955 v19	42,4 ± 4,5
SAFE® 260 HF U8978 v19	43,8 ± 3,4

The animals on the SAFE® 260 HF and the SAFE® 246 HF were both obese with a mean weight respectively of 42,4 ± 4,5 and 43,8 ± 3,4 compared to 28,6 ± 1,6 for the standard SAFE® A03 group (Figure 1). Indeed, we have seen that the significant weight difference appears after 1 week of diet for the SAFE® 260 HF and after 3 weeks for the SAFE® 246 HF but both diets show a great weight gain with similar slope.

We also performed analysis of fat mass of mice (Figure 2). We saw that both high fat diets were significantly different from control but neither the volume of total adipose tissue nor the surface of subcutaneous (scWAT) and intra-abdominal (iaWAT) white adipose tissue were different between SAFE® 260 HF and SAFE® 246 HF.


Figure 2 : Effect of High Fat diets on the increase of adipose tissue seen with tomography experiment. After 12 weeks of diet, mice (n= 10) are put to sleep and then places in a positron emission scanner to obtain 3D reconstitution of the animal and calculate fat and lean mass. (a) Volume of adipose tissue in the whole body of mice. (b) Surface of scWAT (subcutaneous white adipose tissue) and iaWAT (intra-abdominal white adipose tissue) at the L4/L5 vertebrae level. The data is ± SEM and Kruskal-Wallis statistic test is performed between each diet. **p<0.01 ; ***p<0.001 ; ****p<0.0001

At the molecular level, we saw the same conclusion for the adipocytes size (Figure 3). High fat diets have larger adipocytes in scWAT than the control but they're not different between each other.


Figure 3 : Effect of High Fat diets on adipocytes size in subcutaneous adipose tissue. Images of subcutaneous adipose tissue slice (10 µm thickness) stained with Hematoxylin/Eosin

We performed Glucose Tolerance Test to see glycemia regulation after injection of glucose and it seems that the regulation of glucose homeostasis is less altered with 12 weeks of SAFE® 246 HF compared to SAFE® 260 HF because the slope looks like the standard SAFE® A03 group (Figure 4).


Figure 4 : Effect of High Fat diets on glucose homeostasis (a) Glucose tolerance test is performed on mice (n=10) after 12 weeks of diet. Mice are fasted for 16 hours and then injected with 1.5g of glucose/kg. Blood sugar level was measured just before injection and then 15, 30, 45, 60, 90 and 120 minutes after. (b) Basal glycemia is measured just after 16 hours fasting. The data is ± SEM and Kruskal-Wallis statistic test is performed between each diet. ****p<0.0001

SAFE® 260 HF and SAFE® 246 HF diets are effective in inducing obesity in mice at the level of weight gain and adipose tissue. Data not presented here also shown that both diets trigger hyperleptinemia and steatosis pattern in histologic analysis of the liver. However, to induce a type II diabetes-like phenotype, only the SAFE® 260 HF diet is effective because SAFE® 246 HF seems to induce hyperglycemia but no alteration of glucose regulation. In conclusion, SAFE® 260 HF is well used to induce obesity and type II diabetes phenotype and SAFE® 246 HF seems then to be less drastic in inducing phenotype in mice and might be more physiological in comparison to human.

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