### Bedding

## **SAFE** BEDDING

## SAFE<sup>®</sup> aspen fine

PRODUCT DATA SHEET Release date: July 2019

#### Definition

ASPEN WOOD GRANULATE

Made from selected premium raw material (trunk wood)

#### **Product Purpose**

Bedding material

#### Directions for Use

STORAGE

Store in a clean and dry place with constant temperature, protected from sunlight

**SHELF-LIFE** 36 months from date of production

#### Particle Size Distribution

Average particle size/particle range				2 mm - 3 mm
Sieve	> 2.5 mm	> 2 mm	> 1.25 mm	> 0.71 mm
Retained	15 - 45 %	50 - 85 %	90 - 100 %	97 - 100%

### Chemical and Physical Properties

SAFE <sup>®</sup> aspen fine	
Color	yellow
Structure	cubic granulate
Ash (850° C, 4 h)	$\sim$ 1 %
Bulk density	~ 195 g/l
pH-value (10 % suspension)	4.5 - 6.5
100 % natural product	

#### Product Presentation

Standard Packaging: 10 kg autoclavable polypropylene bags.

Big bags and bulk trucks are also available.

Customized solutions and irradiation are possible upon request.

The batch number system provides full traceability.

Analysis certificate of contaminants provided upon request.



SAFE<sup>®</sup> apen fine

#### Certificates









Pictures indicative only

Produced in Germany

# COMPLETE CARE COMPETENCE

WORLDWIDE HEADQUARTERS 73494 Rosenberg (Germany) service@safe-lab.com PRODUCTION SITE DIETS 89290 Augy (France) info@safe-lab.com

### www.safe-lab.com