Bedding



# SAFE<sup>®</sup> performance

PRODUCT DATA SHEET Release date: July 2020

Picture indicative only

#### Definition

CELLULOSE PELLETS Made from selected premium raw material

#### **Product Purpose**

Bedding material

### Directions for Use

**STORAGE** Store in a clean and dry place with constant temperature, protected from sunlight **SHELF-LIFE** 36 months from date of production

# Particle Size Distribution

Average particle size/particle range

6 mm - 15 mm

# Chemical and Physical Properties

SAFE <sup>®</sup> performance	
Color	light grey
Structure	pellets
Bulk density	~ 150 g/l
pH-value (10 % suspension)	5.5 - 7.5
100 % natural product	

Product Presentation

Standard Packaging: 35 l autoclavable paper bags.

Customized solutions and irradiation are possible upon request.

The batch number system provides full traceability.

Analysis certificate of contaminants provided upon request.



SAFE<sup>®</sup> performance

#### Certificates





Produced in USA

PRODUCTION SITE DIETS 89290 Augy (France) info@safe-lab.com

# www.safe-lab.com