Bedding

C S A F E D D I M G

SAFE[®] pure

PRODUCT DATA SHEET Release date: July 2019

Definition

CELLULOSE FLAKES Made from selected premium raw material

Product Purpose

Bedding material

Directions for Use

STORAGE

Store in a clean and dry place with constant temperature, protected from sunlight

SHELF-LIFE 36 months from date of production

Chemical and Physical Properties

| SAFE [®] pure | |
|----------------------------|-----------|
| Color | white |
| Structure | flakes |
| Ash (850° C, 4 h) | max 2.5 % |
| Bulk density | ~ 130 g/l |
| pH-value (10 % suspension) | 5.5 - 7.5 |
| 100 % natural product | |

Product Presentation

Standard Packaging: 7 kg autoclavable polypropylene bags.

Customized solutions and irradiation are possible upon request.

The batch number system provides full traceability.

Analysis certificate of contaminants provided upon request.



SAFE[®] pure

Certificates







Picture indicative only

Produced in Germany



WORLDWIDE HEADQUARTERS 73494 Rosenberg (Germany) service@safe-lab.com PRODUCTION SITE DIETS 89290 Augy (France) info@safe-lab.com

www.safe-lab.com