## Bedding



# SAFE® select fine

PRODUCT DATA SHEET Release date: July 2019

## Definition

#### POPLAR WOOD GRANULATE

Made from selected premium raw material (trunk wood)

## **Product Purpose**

Bedding material

## Directions for Use

#### STORAGE

Store in a clean and dry place with constant temperature, protected from sunlight

#### SHELF-LIFE

36 months from date of production

## Particle Size Distribution

Average particle size/particle range			2.5 mm - 3.5 mm	
Sieve	> 2.5 mm	> 2 mm	> 1.25 mm	> 0.71 mm
Retained	25 - 50 %	60 - 85 %	min. 90 %	min. 97 %



SAFE® select fine

Picture indicative only

## Certificates









## Chemical and Physical Properties

SAFE® select fine	
Color	yellow
Structure	cubic granulate
Ash (850° C, 4 h)	~ 0.5 %
Bulk density	~ 195 g/l
pH-value (10 % suspension)	4.5 - 6.5
100 % natural product	

## **Product Presentation**

Standard Packaging: 10 kg autoclavable polypropylene bags.

Big bags and bulk trucks are also available.

Customized solutions and irradiation are possible upon request.

The batch number system provides full traceability.

Analysis certificate of contaminants provided upon request.



Pictures indicative only



