# Scientific Diets



PRODUCT DATA SHEET

Release date: August 2020

Page 1/2

# SAFE® A30

# Definition

Complete breeding vegetal diet for rats, mice and hamsters.

# **Product Purpose**

Diet for growing and breeding, pregnant and nursing animals. To be used within the context of experimental protocols. Does not contain animal proteins, alfalfa and its byproducts.



### DISTRIBUTION

#### Period

From birth onwards. A transition period to SAFE A40 maintenance diet during weaning is recommended.

#### Method

- Ad libitum or rationed according to experimental protocols.
- Remove from the packaging and place directly in the cage feeder or on the cage floor.
- Keep fresh water always available.

### DAILY CONSUMPTION

Rats 18 to 25 g, mice 3 to 6 g, hamsters 8 to 12 g.

#### **STORAGE**

Store in a clean, dry and cool place, protected from light.

# SHELF-LIFE from the date of production

Paper bag or plastic pouch = 12 months Vacuum packed = 24 months

SAFE® A30

Picture indicative only

# Irradiation

Possible doses: Minimum 10, 25 or 40 kilograys.

# **Product Form**

PELLETS	Mean
Diameter	11.1 mm
Crushing resistance	13.8 kgf/cm <sup>2</sup>
Abrasion resistance	98.9 %
Specific mass	680 g/l
Average pellet weight	2.1 g
Average pellet length	20.2 mm

Also available powdered on demand.

# **Product Presentation**

\*All SAFE® diets are available with different packaging, irradiation and with analytical data on demand. Selected solutions of the most sold items from the SAFE® portfolio.

DIET STANDARD PACKAGING		USUALLY AVAILABLE WITH IRRADIATION DOSE	
SAFE® A30	1 x 10 kg Paper bag		
SAFE® A30 SP*	1 x 10 kg Paper bag in plastic pouch	Min. 10 kGy, Min. 25 kGy	
SAFE® R30*	2 x 5 kg Paper bag, double vacuum packed and boxed	Min. 25 kGy	
SAFE® R30*	10 x 1 kg Double vacuum packed and boxed	Min. 40 kGy	

Produced in France



# Scientific Diets



PRODUCT DATA SHEET

Release date: August 2020

Page 2/2

# SAFE® A30

# Ingredients

Wheat, soybean meal, maize, wheat bran, barley, extruded soybeans, maize gluten, wheat germ, sunflower seed, inactivated brewer's yeast, calcium carbonate, pre-mixture of vitamins, dicalcium phosphate, pre-mixture of minerals, L-lysine, sodium chloride.

# Analysis End Product TOTAL PER KG

### **AMINO ACIDS**

Arginine	12 800 mg	Méthionine	4 100 mg
Cystine	3 800 mg	Tryptophane	2 500 mg
Lysine	13 300 mg	Glycine	12 300 mg


#### **FATTY ACIDS**

Palmitic acid	5 500 mg
Stearic acid	1 600 mg
Oleic acid	14 800 mg
LA	15 000 mg
ALA	1 400 mg

**CENTESIMAL COMPOSITION** 

Cereals	59.5 %
Vegetal Proteins	35.0 %
Vitamins & Minerals	5.1 %
Carbon Hydrates	0.40 %

# NUTRITIONAL COMPOSITION



## **ENERGY CONTENT**

	MJ/kg	kcal/kg	%
ME Pig	13.4	3 208	
ME Atwater	14.2	3 389	
Energy from proteins	3.7	880	26.0
Energy from lipids	2.0	477	14.1
Energy from NFE	8.5	2 032	60.0
More information on energy calculation: www.safe-lab.com			_

**MINERALS END PRODUCT** Calcium 9 500 mg Phosphorus 6 000 mg Sodium 2 500 mg 8 000 mg Potassium Magnesium 1 800 mg Manganese 75 mg Iron 290 mg Copper 20 mg Zinc 68 mg Chlorine 4 800 mg

VITAMINS	END PRODUCT
Vitamin A	14 000 IU
Vitamin D3	2 000 IU
Vitamin E	50 IU
Vitamin K3	5.7 mg
Vitamin B1	8.0 mg
Vitamin B2	13 mg
Vitamin B3	90 mg
Vitamin B5	15 mg
Vitamin B6	3.5 mg
Vitamin B9	0.50 mg
Vitamin B12	0.020 mg
Biotin	0.10 mg
Choline	2 100 mg

For the welfare of animals SAFE® bedding and environmental enrichment such as SAFE® block gnawing logs and SAFE® nesting materials should be available in the cage.

The values of the end products are given as indication only and have no contractual value. They are calculated averages of product analysis results before irradiation and autoclaving. Depending on production conditions, storage and analytical methods variations may occur. An analysis is performed on request.

Produced in France

