## Enrichment



# SAFE® tube small

PRODUCT DATA SHEET Release date: January 2021

#### Definition

CARDBOARD TUNNEL

#### **Product Purpose**

Enrichment material

USE

Nesting	Gnawing	Activity
		Y



#### Directions for Use

STORAGE

Store in a clean and dry place with constant temperature, protected from light

SHELF-LIFE

36 months from date of production

**AUTOCLAVABLE** 



SAFE® tube small

Picture indicative only

### **Product Properties**

SAFE® tube small	
Color	light brown
Structure	tunne
Dimensions	Ø 50 mm·ca 75 mm (H)

#### Certificates







#### **Product Presentation**

Standard Packaging: 300 units per carton

Full pallets are also available.

Customized solutions and irradiation are possible upon request.

The batch number system provides high traceability.

Analysis certificate of contaminants provided every quarter upon request.