

SAFE® compact nesting large

PRODUCT DATA SHEET
Release date: October 2019

Definition

SELECTED SOFTWOOD SHAVINGS, COMPACTED

Product Purpose

Enrichment material

USE

| Nesting | Gnawing | Activity |
|---------|---------|----------|
| x | | x |



Directions for Use

STORAGE

Store in a clean and dry place with constant temperature, protected from light

SHELF-LIFE

36 months from date of production

AUTOCLAVABLE

Product Properties

SAFE® compact nesting large

| | |
|------------|----------------------------|
| Color | light brown |
| Structure | compacted shavings |
| Dimensions | Ø ca. 50 mm; ca. 10 mm (H) |

SAFE® compact nesting large

Picture indicative only

Certificates



Product Presentation

Standard Packaging: 500 units per carton

Full pallets are also available.

Customized solutions and irradiation are possible upon request.

The batch number system provides high traceability.

Analysis certificate of contaminants provided every quarter upon request.